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Numer części: 730022, 730023, 443300

Farnsworth D15 Color Test
- **730022**

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Bibliografia:

1. Birch, J. Diagnosis of Defective Colour Vision. Oxford Medical Publications 1993
2. Farnsworth, D. The Farnsworth Munsell 100-Hue and dichotomous tests for colour vision. Journal of Ophthalmology Society American 33:568, 1943
3. Fukami, K. Evaluation of the Farnsworth-Munsell 100-Hue Test. Japanese Journal of Clinical Ophthalmology 30:27-31, 1976
4. Greenstein V, Sarter B, Noble K, and Carr R. Investigative Ophthalmology & Visual Science Vol. 31, 1008-1014. Hue discrimination and S cone pathway sensitivity in early diabetic retinopathy.
5. Hahn C. Evaluation of Hahn Double 15 Hue Test. Poster Session 10th Japan-Korea Joint Meeting of Ophthalmology 9/21-23/2000
6. Helve, J. A comparative study of several diagnostic tests of colour vision used for measuring types and degrees of congenital red-green defects. Acta Ophthalmologica Supplement. 115:18, 1972

7. Hyvärinen L. Quantitative Color Vision Test V-16 Manual.
8. Smith VC, Pokorny J. Large-field trichromacy in protanopes and deuteranopes. Journal of Opt Soc Am 1977; 67:213-220.
9. Tasman W, Jaeger E.A. Duane’s Clinical Ophthalmology, Vol 3, Chapter 6, Lippincott Williams & Wilkins, 2000
10. Verriest G, van Laethem J, Uvijel A. A new assessment of the normal ranges of the Farnsworth-Munsell 100 hue test scores. American Journal of Ophthalmology 1982; 93:635-642
11. Vingrys AJ and King-Smith PE. Investigative Ophthalmology & Visual Science Vol. 29, 50-53. ‘A quantitative scoring technique for panel tests of color vision’.
12. McIntyre, Donald. Colour Blindness - Causes and Effects. Dalton Publishing, 2002

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1. Bernstein HN. Ophthalmic considerations and testing patients receiving long-term antimalarial therapy. AM J Med 1983;75(18):25-34
2. Muirhead KD. The use of chloroquine and D-penicillamine in the treatment of rheumatoid arthritis. Med J Aus 1986;144(1):32-7.
3. Cullen AP, Chou BR. Keratopathy with low dose chloroquine therapy. J Am Optom Assoc 1986;5(5):368-72
4. Meischer PA. Treatment of systemic lupus erythematosus. Springer Sem Immunopathol 1986;9:271-82
5. Bartlett JD, Jaanus SD. Ocular effects of systemic drugs. In: Bartlett JD, Jaanus SD, eds. Clinical Ocular Pharmacology, 2nd ed. Boston: Butterworths, 1989: 901-42
6. Finbloom DS, Silver K, Newsome DA, Gunkel R. Comparison of hydroxychloroquine and chloroquine use and the development of retinal toxicity. J Rheumatol 1985;12(4):692-4
7. Johnson MW, Vine AK. Hydroxychloroquine therapy in massive total doses without retinal toxicity. 1987;140:139-44
8. The Medical Letter on Drugs and Therapeutics. 1987;29(734):21-4
9. Rynes RL. Ophthalmologic safety of long-term hydroxychloroquine sulfate treatment. Am J Med 1983;75(18):35-8
10. Kastrup EK, et al., eds. Facts and Comparisons. St Louis. Lippincott, 1986:253e-g.
11. Bowman KJ. A method for quantitative scoring of the Farnsworth Panel D-15. Acta Ophthalmologica (Copenh) 1982;60:907-16. Richmond Products Inc.
12. Bowman KJ, Collins MJ, Henry CJ. The effect of age on performance of the panel D-15 and de-saturated D-15: a quantitative evaluation. In: Verriest G, ed. Colour Vision Deficiencies, VII. The Hague: Dr W Junk Publishers, 1984:227-31
13. Cyert L. Eye and Vision Conditions in the American Indian. Pueblo Publishing Goss and Edmonson 1990:137-147

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NZHTA Report 7 - New Zealand Health Technology Assessment (NZHTA) The Clearing House for Health Outcomes and Health Technology Assessment Department of Public Health and General Practice Christchurch School of Medicine Christchurch, N.Z. Colour vision screening A <http://nzhta.chmeds.ac.nz/colour.htm#screening> <http://orlab.optom.unsw.edu.au/ICVS/Daltoniana.April98.html>

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Numer części: 730022, 730023, 443300

1155 Jansen Farm Drive
Elgin, IL 60123 USA

Toll Free: 800-362-3860 Phone: 847-841-1145
Toll Free Fax: 888-362-2576 Fax: 847-841-1149

www.good-lite.com


Good-Lite B.V.
Kantstraat 19, NL-5076NP Haaren
P.O. Box 10110, NL-5060 GA Oisterwijk
The Netherlands
EU@medical-risk.com
+31 (0) 411.623.791
SRN: NL-AR-000002695

QNET Ltd.
Livingstone House, 309 Harrow Road
Wembley, Middlesex HA9 6BD (GB)
UK@medical-risk.com
MHRA Account No: 0000011991

QNET CH-REP GmbH
Im Büel 15
8750 Glarus (Switzerland)
CH-REP@medical-risk.com
UID Number: CHE-313.448.831
Ph.: +41 79 610 0870

July 2024

Date of issue: 18.08.2025

Revision: 0.01

MD Farnsworth D15 Color Test

REF 730022

MD Lanthony Desaturated Color Test

REF 730023

MD Lanthony D-15 Color Test(Magnetic) - Replacement Color Chip

REF 443300