



LEA SYMBOLS Near Screener and Distance Screener Book

Part Number: 252000 & 259900

- A. Near Vision Screener with 16 inch/40 cm Cord (#252000)
- B. Distance Vision Screener Book, 10 feet/3 meter (#259900)



The most effective vision screening for detecting lazy eye uses near and distance test. When testing children younger than 4 years, testing is a lot easier if children have had an opportunity to play with the test symbols before the test situation. This can be arranged by including the test symbols in the information sent or given to parents as an invitation to participate in the screening examination. Parents can play the *Lea Symbols* game at the *Lea-Test* homepage. Just prior to the testing children can play with the *Lea Puzzle* that effectively decreases fears and introduces the test situation as a play situation.

Near vision is functionally more important than distance vision in the life of a young child. The child is also more accustomed in using vision at near than at greater distances. Therefore introduction of the test at near familiarizes the child with the test situation. Young children are easier to test at near than at distance of 3 meters (10 feet). In the rare case of myopia you will find that the child has useful vision at near and parents (and you) will not be alarmed when the child does not see well during the distance test.

The target population of amblyopia screening, like any screening, are symptom free children. If a child has symptoms of any kind (strabismus, squeezing eyes during testing, head tilt or turn, itchy or red eyes etc.), (s)he is referred for treatment independent of visual acuity values.

The *Lea Symbols Amblyopia Screeners* look different than the other *Lea Symbols* test. The grey background was chosen for two reasons: 1. in developing countries dust will not show on the test during the day, in the evening the tests can be washed, 2. since there is less white surface, the test does not dazzle photophobic children.

The test lines are wider apart in the near vision screener card than in the standard test (#250800). It makes it easier for both the tester and the child to know on which line to read. Therefore testing becomes more relaxed and is faster. In the distance tests there is only one line visible on each page and there are two pages of each size. This way if the child makes several errors on a certain line but immediately corrects some of them, the tester can ask the child to read the other line of the same size.

Instructions

- Establish a method of communication such as naming (signing) or pointing (matching). Decide with the child which names will be used to identify the symbols. When needed, train with the *LEA Puzzle* (#251600) or the *Response Key Card* (#251700).
- Start with binocular testing at near. Point to each of the four symbols (circle, house, apple, square) on the top line of the near card; observe the baseline responses for comprehension, speed and accuracy.
- Ask the child to identify only the first symbol on each line.
- Repeat this procedure for each or every second line (moving quickly down the near card to avoid tiring the child; in distance testing move to the next size, flip two pages) until the child hesitates or misidentifies a symbol.
- Move back up one line and ask the child to identify all the symbols on that line.
- If the child identifies all symbols correctly go to the next line down and ask the child to identify all the symbols on that line. The child may have focused the image more carefully and therefore can read symbols that were not seen a moment earlier.
- If the child skips a symbol, ask the child to try again while briefly pointing to that symbol. Do not leave your finger or a wooden stick (never a pen or pencil) at the symbol because it gives a visual reference and may improve fixation of an amblyopic eye.

If the child does not function well in the test situation at near, let him/her train by using the training test that you can print from the *Section Info for Parents* found on the *Lea-Test* website (<http://www.lea-test.fi>). If after 2-3 weeks of training a three year old “does not co-operate”, the child needs to be referred. The reason for poor co-operation may be that the child does not see the symbols! Therefore before a diagnosis of delayed intellectual development is made, it is mandatory to know that sensory functions are normal. This includes not only visual acuity but also other visual cortical functions and hearing.



LEA SYMBOLS Near Screener and Distance Screener Book

Part Number: 252000 & 259900

Testing at distance of 3 meters (10 feet)

After near testing introduce the distance chart to the child by saying, "Let's look at the same pictures a little further away." Move the chart gradually back to 3 meters (10 feet), while watching the child for signs of inattention. If the child loses interest move closer to 2 meters (80 inches), or one meter (40 inches). Always test well within the visual sphere of the child. Older children may be switched directly from a near vision test to the 3 meter chart.

After binocular testing, proceed with testing each eye separately. Use two pairs of plano glasses or the special screening frame for occlusion of the child's eyes, one at a time. This is the least disturbing type of occluder for children.

For monocular testing, follow the same procedure as for binocular testing.

If the near vision card is held at 40cm (16 in) and the distance visual acuity test at 3 meters (10 feet) the visual acuity value is found in the margin adjacent to that line. The visual acuity is recorded as the last line on which at least 3 of the 5 symbols are read correctly.

Always test until the threshold line so that you can measure the difference in visual acuity between the eyes.

The tests are supplied with a response key. For convenience, the near vision card has a response key line at the bottom. However, naming is a much faster way of responding and you can hear when the child becomes hesitant or starts to lose interest.

Monocular near vision testing

As a part of vision screening, monocular testing at near is often possible much earlier than testing at distance. In testing young children this may be the only way of measuring monocular values.

Monocular testing at near is of special interest at the age of 6-7 years: If visual acuity at distance has become less than it was at age 4 or if there is difference between the two eyes, near vision measurement may give the following diagnosis: If the near vision values are symmetric and better than the values at distance, the change in distance vision cannot be caused by anything else but mild myopia, which does not need to be corrected. The child does not need to be referred. This leads to a decrease in the expenses of vision screening and simultaneously leads to an improvement in the quality of screening.

If there is decrease in both the distance and the near vision acuity, the most common reason for that is a change in refractive error but it may be caused by a disease. Therefore the child needs to be referred.

Pass/fail criteria

Pass/fail limits are used in a number of amblyopia screenings. It decreases the sensitivity of the test.

If the child has visual acuity 0.8 (20/25, 6/9) in one eye and 0.5(20/40, 6/12) in the other eye, the difference in visual acuity values between the eyes will not be detected if a pass/fail criteria of 0.5 (20/40, 6/12) is used.

At the age of four years, most western countries use a pass/fail criteria of 0.5 (20/40, 6/12) for binocular visual acuity to detect large refractive errors and previously undetected cases of binocular visual impairment that might affect the child's general development.

On the other hand, a pass/fail limit of 0.8 (20/25, 6/9) at the age of six years does not seem to pick many children who have refractive errors that should be corrected.

In school age the passing limit should not be one single visual acuity value but based on the vision needs of a particular child: If the child is in the shortest or the middle third of the class and if (s)he can sit close to the podium, binocular visual acuity of 0.2 - 0.3 (20/100 - 20/60, 6/30 - 6/18) may cause no difficulty in the class room. On the other hand if the child is tall and is required to sit far from the blackboard, visual acuity of 0.8 (20/25, 6/9) may be required for normal working in the classroom. If children are referred at a time when they have only mild myopia, the referral expenses are high. In many places children are prescribed these weak glasses that they do not need. This is a sizeable extra cost to the families.

It should be clearly remembered that vision screening after the age of 8-9 years is arranged to find those children who may need glasses in order to see well in the classroom. We are not looking for amblyopia or diseases and therefore the deciding fact is whether the child sees well enough in the classroom or not.

There is one more question that needs to be stressed: age normal visual acuity values in vision screening do not guarantee that vision is normal. Brain damage related vision loss (CVI), retinal degeneration or optic nerve lesions may not affect visual acuity at high contrast but may have caused loss of vision at low contrast or visual field defects or disabling losses of visual perception. The symptoms of brain damage and eye or pathway disorder related changes should be better known because today it may take years before a child is referred to intervention and receives services of special education.

252000 - Answer Key

Testing distance of 10ft/3m

Distance Equivalents

Foot	Meter	Front of Card	Letter Size	Decimal
20/400	6/120	🏠 □ ○ □ ♥	8.0M	.05
20/320	6/95	○ 🏠 ♥ ○ □	6.3M	.063
20/250	6/75	🏠 □ 🏠 ♥ ○	5.0M	.08
20/200	6/60	□ 🏠 ♥ ○ ♥	4.0M	.10
20/160	6/48	🏠 □ ○ ♥ □	3.2M	.12
20/125	6/38	○ ♥ □ 🏠 ○	2.5M	.16
20/100	6/30	♥ ○ 🏠 □ 🏠	2.0M	.20

Back of Card

20/100	6/30	♥ □ ○ □ 🏠	🏠 □ ○ 🏠 ♥	○ □ ○ 🏠 ♥	2.0M	.20
20/80	6/24	□ ○ ♥ 🏠 ○	□ ♥ ○ 🏠 □	□ ○ 🏠 ♥ □	1.6M	.25
20/63	6/19	○ ♥ 🏠 □ ○	♥ 🏠 ♥ □ ○	🏠 □ ♥ 🏠 ○	1.25M	.32
20/50	6/15	♥ ○ ♥ 🏠 □	🏠 ○ □ ♥ 🏠	○ □ 🏠 ○ ♥	1.0M	.40
20/40	6/12	□ ♥ ○ □ 🏠	○ □ ○ 🏠 ♥	♥ 🏠 □ ○ 🏠	.80M	.50
20/32	6/9.5	○ 🏠 □ ♥ ○	□ ○ □ ♥ 🏠	□ ♥ 🏠 □ ○	.63M	.63
20/25	6/7.5	🏠 □ 🏠 ○ ♥	○ □ ♥ 🏠 ○	🏠 ♥ ○ ♥ □	.50M	.80
20/20	6/6	□ 🏠 ○ ♥ □	🏠 ♥ □ ○ ♥	○ □ ♥ □ 🏠	.40M	1.0
20/16	6/4.8	○ □ ♥ 🏠 ♥	♥ 🏠 ○ □ ○	🏠 ○ □ 🏠 ♥	.32M	1.25